



ECL MP 121/151 MULTI-PARAMETER MONITOR

EXQUISITE AND CONVENIENCE

12.1 or 15 inch high brightness long-life LED display, optional touch screen

Lightweight and quick-hook design, simple fix on poles less than 30 mm in diameter

Iconify key operation, centralized menu make it easier and faster to use



FLEXIBLE EXTENSION INTERFACE, FREE SELECTION

One-key switch and indenting key menu save more display space

Optional functions: IBP, EtCO₂, NELLCOR or MASIMO SPO₂, SUNTECH NIBP, REC, touch, HDMI, WIFI, Central Station software, Scanning gun, Card reader and information input module ect

PARAMETERS

Model	EGL MP 121	EGL MP 151
Screen size	12 inch	15 inch
ECG	Standard three or five lead cable	
HR	Adult 15~300 bpm, Pediatric 15~350 bpm	
RESP	Range: adult 2~120 BrPM, pediatric: 2~150 BrPM; Accuracy: ± 2 BrPM	
TEMP	Range 0~50°C Accuracy: $\pm 0.1^\circ\text{C}$	
PR	Range 20~300 bpm, Accuracy: ± 3 bpm	
SPO2	Range 0~100%, Accuracy: 70~100% ± 2	
NIBP	Operating mode: Manual / Automatic / STAT Range 10~270 mmHg, Accuracy: ± 5 mmHg	
Interface	Ac power connector, USB, RJ45, HDMI etc	
Standard configuration	ECG cable, ECG electrode, NIBP cuff, NIBP tube, SpO2 probe, Temperature probe each one	
Optional features	IBP, EtCO2, NELLCOR or MASIMO SPO2, SUNTECH NIBP, REC, TOUCH, WIFI, Scanning gun, Card reader and information input module etc	

ACCESSORIES (OPTIONAL)

			
SPO2 Probe	NIBP Tube	IBP Cable	Temperature Probe
			
ECG cable	NIBP Cuff	EtCO2	EtCO2

HILD LIFESCIENCES PRIVATE LIMITED

Eagle Tower, 2nd Floor, 27, Greams Road, Thousand Lights, Chennai - 600 006, Tamil Nadu, INDIA
 Ph: +91(44) 4324 1900 - 09 | email: info@hildmedical.com | website: www.hildmedical.com